

MC MEDICREATIONS

The MediShape



**Safe, Proven, Aggravation-Free
Body Contouring.**



No Rollers Or Consumables. One Cable.

A MediShape multipolar RF device can be a potent solution to cellulite, loose skin, and that stubborn fat around the waist.

We designed it with two major goals in mind. The first—to maximize effectiveness while minimizing risk. The second?

Zero annoyances.

That means we skipped out on the unwieldy cameras. We resisted adding in awkward straps. The handpieces are light enough that a 30-minute session won't actually give you a workout.

What's left is intuitively simple (but powerful) bipolar RF energy treatment that leaves both you and your patients happy with the results.

America's Leading Clinics Use Medicreation Devices.

Before

After



Before

After



The MediShape:

Product	Multi-Functional Shape Device
Vacuum	Max Vacuum: 350mm Hg
RF Frequency	1.2M Hz Bipolar, Multi Polar, Max 50W
LED Light Therapy	Red LED (630nm, Blue LED (430nm)
Input Power	AC 110V
Vacuum Modes	Continuous, Pulse 1-7, Super Pulse

Sick of the Status Quo? So Are We.

The Medicreations team has nearly

50

years of combined experience immersed in some of the world's leading aesthetic device brands.



We've learned how things are "supposed to go" from the giants of the industry. **It's terrible.**



Devices built around... Expensive, single-use consumables. Manipulative sales tactics. Absurd pricing structures designed to gouge every last cent from care providers.



So... We got together some pretty smart scientists, engineers, and business people to do things differently.

Our revolutionary idea:

Make a living *with* the people who use our devices, not *from* them.

Quality Devices + Transparent Process + Fair Pricing = Our Promise

Aesthetic care professionals should have access to reliable, high-performance devices at prices that actually make sense.

That's why we design our machines to utilize *as few consumables* as possible, even when competitors don't. Why we have a *industry-leading two-year guarantee* for all parts and damage. Why even though we'd put our devices against nearly any other brand's in terms of performance, *we don't charge a premium*.

As for our process, our representatives will never sell you a device that you do not need. If it won't improve your business and make your life richer, it doesn't go out the door. Period.

“The way it's always been” isn't the right way. We want to do better.

Always Built In Our FDA-Compliant Facility Right Here In Las Vegas, NV.

Medicreations devices are some of the only aesthetic care devices fully manufactured in the United States, with at least

61%


of all components sourced from U.S. suppliers.

Our FDA-compliant facility features a clean room and full-service center, staffed by our excellent team of professionals. Every device is meticulously built and subject to strenuous quality control checks before shipping out.


We film videos of our process and offer guided tours of our facility on request. We are proud of how we build our devices and what goes into them.

If you're ever down in Vegas, let us know and we'll show you around!

How Does A MediShape Treatment Burn Fat And Tighten Skin?



RF energy reaches deep into the skin.




Energy flips water molecules to create friction.




Heats up fat, stimulates fibroblasts, and increases circulation.

Smooths Fat



Heat liquidizes fat, making it more malleable.

Tightens Skin



Stimulated fibroblasts boost collagen and elastin production.

Removes Fat



Increased blood flow induces redistribution and expulsion of fat.

Made For The Working Practitioner— Ergonomically Superior. Faster.

Better Ergonomics.

A single cable means no tangling, no tripping. The handpieces are light and fit naturally in the hand during usage. They were meticulously designed to be the perfect size. Small enough to be nimble. Large enough to pack a punch.

No Gimmicky Addons.

Say goodbye to the roller—few will miss it. Its replacement is a set of multipolar handpieces that are safer and easier to use. The MediShape includes vacuum suction, ultrasound cavitation, and LED light therapy to supplement and enhance treatment.

Faster Treatments.

Multipolar RF technology means more contact with the skin and more energy utilized per second. Therefore, the device heats the affected tissue to the required 40°C - 45°C far quicker, leading to faster and more effective sessions.

No Consumables.

You don't need to buy any kind of special disposables from us to get top-notch results. Simply use your preferred gel or mineral oil and get on with your day.

Easier Cleaning.

No rollers means less chance of gunk getting stuck between the cracks. A few minutes it all it takes between sessions to have a MediShape looking picture-perfect.

Built-In Patient Comfort And Safety. Less Arcing And Pinching. Maximum Effect.

The MediShape's bipolar RF technology hands you absolute control over how deep the energy penetrates, allowing you to customize treatment for each patient and keep discomfort minimal (or non-existent).

A 350mm Hg vacuum with multiple settings (continuous, pulse 1-7, and super-pulse) helps ensure that, during use, the RF energy probes stay in full contact with the skin. Multiple points of contact boosts the effectiveness of the treatment while also lowering the chance of arcing energy.

There are three fully detachable RF handpieces available. Small, medium, and large. This guarantees you always have the perfect size to maximize speed of treatment while ensuring comfort for both you and your patient.

An ultrasound cavitation handpiece is also included to tighten skin while stimulating collagen and elastin, providing a truly comprehensive body contouring experience.

This is safe, effective fat burning and body contouring packaged into a streamlined device **you'll actually like using.** Without the complications.



Every Medicreations Device Is Backed By An Industry-Leading 2-Year “We’ll Fix It” Warranty.

In the rare event your devices malfunction or break due to a technical defect—**we’ll take care of it.**

The benefit of having our manufacturing facility in the U.S. is that we always have access to the very same engineers and technicians who designed and built your machine in the first place. We know exactly how to get your devices working perfectly again. And if we don’t? We’ll replace it at no charge. **Simple.**

We stand behind our products for a full two years—double that of the industry standard. Nor do we require you to pay tens of thousands of dollars for that warranty.

If you buy any Medicreations device, you just won’t have to worry about it for a full 730 days, with options past that.



The Whole Is Greater Than The Sum Of Its Parts, And **Definitely** One Part Alone.

To most med spa operators, a single device can seem like an incredible investment. Here's the calculus.

If a machine works out well, the buyer might earn a return in 2-3 years. At which point they might buy a different one. Rinse and repeat, and they could one day provide the full range of experiences their patients are looking for. Theoretically.

Let's skip that.

Medicreations offers a whole suite of products to deliver comprehensive aesthetic treatments for patients, all at the cost of a single device from other manufacturers. You can provide hair removal alongside skin tightening, body sculpting with muscle toning, or any other synergistic combination. No waiting for years—you could own them all tomorrow.

The difference between 1 device & several is incalculable.

It's a single instrument against an orchestra, a solo competing with a symphony. Only when all the pieces meld together will you find yourself producing that cohesive experience that is "just so."

We'll Teach You Everything About Your Device (And More Besides).

We *love* helping people learn the science and use behind our products. Seriously. So much so we actually built a little academy here in Vegas to do it better.

You'll receive a comprehensive onboarding process when you get any of our lineup at no additional cost. Everything you need to know about usage and care of your devices will be included. That's expected, a given. But we offer more.

Regular Webinars With Top Doctors On The Science

In-Person Seminars With Use-Cases

Live Demonstrations

Q&A Sessions

Marketing Help

Whatever you want to learn, we'll provide.

See If A Medicreations Device Is Perfect For You

Lets talk. We'll make sure you understand the science behind the device you're interested in, who it's ideal for, and how to get the most from it. And if during the process, you find out you're not so sure it's for you?

That's OK.

We only want to sell devices that we sincerely believe will be a benefit to your practice, help you hit the goals you want to hit, and blow your ROI out of the water.

Something that checks every box for you in particular, with NO exceptions.

Maybe that something will be a Medicreations device.

Contact:

Private Practice Solutions LLC

www.ppslasers.com

E: info@ppslasers.com | T: (323) 389-6028

Disclaimer: No representation or warranty will be implied from any description of or claims regarding the equipment or its effectiveness or the ability to achieve any particular clinical results, whether written or oral, as contained in this brochure or similar statements made or furnished to buyer by any person including Medicreations personnel. The use of the equipment requires the expertise of sound professional judgment and the treatments results may be based on operator skill and experience, patient/client suitability, patient/client response to treatment and other factors beyond the control of Medicreations. Operation of the equipment must be in accordance with all applicable laws and regulations and medical and treatments guidelines and for ensuring that each operator of the equipment is adequately trained and qualified to use and operate the equipment safely and properly and to perform procedures in accordance with such law, regulations and guidelines. Medicreations makes no representation or warranties regarding federal, state or local laws or regulations, or medical or treatments guidelines that may apply to the use and operation of the equipment. Published on March 13th, 2020.