

The MediSpa



The Ultimate All-In-One Device For Estheticians.



Every Medical Spa Needs The True Facial Cleanse™ In Their Toolkit.

The MediSpa packs practically everything you can do for someone's face into a single device, without reaching for a laser or needle.

It exfoliates (twice over). It infuses serums. It tightens skin. It clears pores. And it absolutely hypercharges collagen and elastin production.

It lets you offer the ideal “beginner's package” for a new customer with immediate results. The True Facial Cleanse™ has no downtime. No pain.

The MediSpa does more jobs than practically any other Esthetic device and does them perfectly.

Every treatment solution on the MediSpa is backed by proven science and has been carefully designed for top-level efficacy and comfort—without compromises. Each handpiece works with the next to maximize session outcomes.

If you're ready to deliver some of the most popular and effective treatments out there without the inconvenience (and cost) of half a dozen pieces of equipment, keep reading to see how a MediSpa can help.

America's Leading Clinics Use Medicreation Devices.

Before

After



"I love the MediSpa, it has a lot of different modalities. Also, it's great for add-on services. Someone may come in for something quick and basic, but to add-on something of this caliber makes such a huge difference. It's definitely the machine of the future."

**Bri, Esthiology Educator
Aveda Insitute**

The MediSpa:

Hydra-Dermabrasion	Yes
Quantity of Serum Bottles	3
Device cleaning material	1
Ultrasound Handpiece	Yes
Bi-Polar RF Handpiece	Yes, 6 Contact Points
Spray Handpiece (i.e. Oxygen)	Yes
Electro ION, Cold & Hot	Yes
Automatic clean system	Yes

Sick of the Status Quo? So Are We.

The Medicreations team has nearly

50

years of combined experience immersed in some of the world's leading aesthetic device brands.



We've learned how things are "supposed to go" from the giants of the industry. **It's terrible.**



Devices built around... Expensive, single-use consumables. Manipulative sales tactics. Absurd pricing structures designed to gouge every last cent from care providers.



So... We got together some pretty smart scientists, engineers, and business people to do things differently.

Our revolutionary idea:

Make a living *with* the people who use our devices, not *from* them.

Quality Devices + Transparent Process + Fair Pricing = Our Promise

Aesthetic care professionals should have access to reliable, high-performance devices at prices that actually make sense.

That's why we design our machines to utilize *as few consumables* as possible, even when competitors don't. Why we have a *industry-leading two-year guarantee* for all parts and damage. Why even though we'd put our devices against nearly any other brand's in terms of performance, *we don't charge a premium*.

As for our process, our representatives will never sell you a device that you do not need. If it won't improve your business and make your life richer, it doesn't go out the door. Period.

“The way it's always been” isn't the right way. We want to do better.

Always Built In Our FDA-Compliant Facility Right Here In Las Vegas, NV.

Medicreations devices are some of the only aesthetic care devices fully manufactured in the United States, with at least

61%

of all components sourced from U.S. suppliers.

Our FDA-compliant facility features a clean room and full-service center, staffed by our excellent team of professionals. Every device is meticulously built and subject to strenuous quality control checks before shipping out.

We film videos of our process and offer guided tours of our facility on request. We are proud of how we build our devices and what goes into them.

If you're ever down in Vegas, let us know and we'll show you around!

7 Incredible Treatments. 1 Device.

Hot/Cold:

On HOT, the handpiece heats up to 42°C, the ideal (but comfortable) temperature for widening pores for deeper cleaning. On COLD, it cools to 5°C to shrink pores, reduce inflammation, and speed recovery.

Enriched-Oxygen Capsule:

Deliver top-tier exfoliation without the risks and “sandpaper” feel of microdermabrasion. Rubbed on the skin, a gel applied prior to treatment reacts to the capsule to create a CO₂-rich foam. This process scrubs the surface and kills bacteria while also boosting blood circulation, cell turnover, and product absorption.

AquaPeel:

The Aqua Peel provides a second layer of exfoliation and comes with three different chemical solutions that respectively smooth, clean, and brighten the skin. It is an unrivaled extraction tool, with its potent vacuum suctioning out oil, dirt, and dead skin cells, while reaching deeper to extract blackheads and cleanse pores.





Electro-ION Microcurrent:

The MediSpa's Electro-ION plumps skin and gives facial muscles a "workout" through direct low-level electrical energy stimulation. This process catalyzes ATP, which speeds tissue rejuvenation while also helping contour the face.



Ultrasound:

An ideal method for tightening the skin around the brow, chin, and neck, it also helps solutions penetrate deep beyond the surface. But its most powerful benefit is subtle. Ultrasound energy hits below the dermis to catalyze collagen and elastin production, enhancing long-term skin health and feel.



Bi-Polar RF (Radio Frequency):

Bi-Polar RF is an extremely controlled method of targeting specific areas on the face for skin rejuvenation. The RF energy heats up the collagen in the target area to around 40-45°C, which triggers the body's natural regenerative processes and energizes production. This smooths out wrinkles, removes lines, and produces a more youthful appearance.



Oxygen Spraying:

Excellent for providing that extra boost of moisture at the end of a major treatment and ensuring the skin has the necessary building blocks to regenerate quickly. Also gives that sought-after "celebrity glow."

Each Handpiece Stands On Its Own. But When Used Together...


We've put meticulous care into ensuring that every handpiece can provide results patients can see and feel when used on their own (or with another device you own).

That said? Every piece has been built to work with the next.

The heating handpiece widens pores to maximize the effectiveness of the Enriched-Oxygen Tablet exfoliation. That, in turn, preps the skin for the Aqua Peel extraction. Microcurrent, Ultrasound, and RF Energy handpieces synergize to turn up the dial on collagen and elastin production. Three different skin tightening treatments are layered over one another for maximum effect.

That's the True Facial Cleanse™.

The MediSpa allows you to offer your patients the ultimate non-invasive facial rejuvenation package.



Every Medicreations Device Is Backed By An Industry-Leading 2-Year “We’ll Fix It” Warranty.

In the rare event your devices malfunction or break due to a technical defect—we’ll **take care of it.**

The benefit of having our manufacturing facility in the U.S. is that we always have access to the very same engineers and technicians who designed and built your machine in the first place. We know exactly how to get your devices working perfectly again. And if we don’t? We’ll replace it at no charge. **Simple.**

We stand behind our products for a full two years—double that of the industry standard. Nor do we require you to pay tens of thousands of dollars for that warranty.

If you buy any Medicreations device, you just won’t have to worry about it for a full 730 days, with options past that.



The Whole Is Greater Than The Sum Of Its Parts, And Definitely One Part Alone.

To most med spa operators, a single device can seem like an incredible investment. Here's the calculus.

If a machine works out well, the buyer might earn a return in 2-3 years. At which point they might buy a different one. Rinse and repeat, and they could one day provide the full range of experiences their patients are looking for. Theoretically.

Let's skip that.

Medicreations offers a whole suite of products to deliver comprehensive aesthetic treatments for patients, all at the cost of a single device from other manufacturers. You can provide hair removal alongside skin tightening, body sculpting with muscle toning, or any other synergistic combination. No waiting for years—you could own them all tomorrow.

The difference between 1 device & several is incalculable.

It's a single instrument against an orchestra, a solo competing with a symphony. Only when all the pieces meld together will you find yourself producing that cohesive experience that is "just so."

We'll Teach You Everything About Your Device (And More Besides).

We *love* helping people learn the science and use behind our products. Seriously. So much so we actually built a little academy here in Vegas to do it better.

You'll receive a comprehensive onboarding process when you get any of our lineup at no additional cost. Everything you need to know about usage and care of your devices will be included. That's expected, a given. But we offer more.

Regular Webinars With Top Doctors On The Science

In-Person Seminars With Use-Cases

Live Demonstrations

Q&A Sessions

Marketing Help

Whatever you want to learn, we'll provide.

See If A Medicreations Device Is Perfect For You

Lets talk. We'll make sure you understand the science behind the device you're interested in, who it's ideal for, and how to get the most from it. And if during the process, you find out you're not so sure it's for you?

That's OK.

We only want to sell devices that we sincerely believe will be a benefit to your practice, help you hit the goals you want to hit, and blow your ROI out of the water.

Something that checks every box for you in particular, with NO exceptions.

Maybe that something will be a Medicreations device.

Contact:

Private Practice Solutions LLC

www.ppslasers.com

E: info@ppslasers.com | T: (323) 389-6028

Disclaimer: No representation or warranty will be implied from any description of or claims regarding the equipment or its effectiveness or the ability to achieve any particular clinical results, whether written or oral, as contained in this brochure or similar statements made or furnished to buyer by any person including Medicreations personnel. The use of the equipment requires the expertise of sound professional judgment and the treatments results may be based on operator skill and experience, patient/client suitability, patient/client response to treatment and other factors beyond the control of Medicreations. Operation of the equipment must be in accordance with all applicable laws and regulations and medical and treatments guidelines and for ensuring that each operator of the equipment is adequately trained and qualified to use and operate the equipment safely and properly and to perform procedures in accordance with such law, regulations and guidelines. Medicreations makes no representation or warranties regarding federal, state or local laws or regulations, or medical or treatments guidelines that may apply to the use and operation of the equipment. Published on March 19th, 2020.