

MC MEDICREATIONS

The MediYag



**Multi-Wavelength Q-Switched
Laser Treatment—Perfected.**



The **Gold** Standard.

Why are Q-Switched YAGs the go-to choice for tattoo removal and hyperpigmentation issues?

Mainly because they have the ability to output extraordinary amounts of energy in truly miniscule timeframes.

Fractions of a second, in fact.


That supremely high-intensity, short duration method triggers a unique photoacoustic effect that helps treat conditions other lasers are powerless against—with minimal chance of long-term damage to the patient.

This makes them ideal not only for tattoo removal, but many kinds of pigmentation conditions, skin rejuvenation, vascular and pigmented lesions, pore reduction, and wrinkles.

Q-Switched YAGs can specifically “target” and break up specific characteristics in a certain area. Like, for instance, ink pigment in someone’s skin, which is then absorbed and removed by the body’s natural processes.

The MediYag is the result of decades of continuous refinement of laser technology. It has top-of-the-line power and speed, along with a streamlined build with a focus on getting results.

America's Leading Clinics Use Medicreation Devices.



The MediYag:

Laser Wavelength (nm)	1064nm/532nm
Single Pulse	1064nm: 100mj-2600mj
Double Pulse	532nm: 100mj-1300mj
Laser Mode	Single/Double Pulse
Pulse Repetition Rate (Hz)	1Hz-10Hz
Spot Diameter (mm)	2mm-10mm adjustable
Power Supply)	1500W

Sick of the Status Quo? So Are We.

The Medicreations team has nearly

50

years of combined experience immersed in some of the world's leading aesthetic device brands.



We've learned how things are "supposed to go" from the giants of the industry. **It's terrible.**



Devices built around... Expensive, single-use consumables. Manipulative sales tactics. Absurd pricing structures designed to gouge every last cent from care providers.



So... We got together some pretty smart scientists, engineers, and business people to do things differently.

Our revolutionary idea:

Make a living *with* the people who use our devices, not *from* them.

Quality Devices + Transparent Process + Fair Pricing = Our Promise

Aesthetic care professionals should have access to reliable, high-performance devices at prices that actually make sense.

That's why we design our machines to utilize ***as few consumables*** as possible, even when competitors don't. Why we have a ***industry-leading two-year guarantee*** for all parts and damage. Why even though we'd put our devices against nearly any other brand's in terms of performance, ***we don't charge a premium.***

As for our process, our representatives will never sell you a device that you do not need. If it won't improve your business and make your life richer, it doesn't go out the door. Period.

“The way it's always been” isn't the right way. We want to do better.

Always Built In Our FDA-Compliant Facility Right Here In Las Vegas, NV.

Medicreations devices are some of the only aesthetic care devices fully manufactured in the United States, with at least

61%

of all components sourced from U.S. suppliers.

Our FDA-compliant facility features a clean room and full-service center, staffed by our excellent team of professionals. Every device is meticulously built and subject to strenuous quality control checks before shipping out.

We film videos of our process and offer guided tours of our facility on request. We are proud of how we build our devices and what goes into them.

If you're ever down in Vegas, let us know and we'll show you around!

5 Nanosecond Pulses. 1300 mJ.

The MediYag laser fires at a 5 nanosecond pulse duration—up to 10hz.

Long enough to bring the full power of its energy upon the target. Short enough to minimize side effects for the patient, regardless of their skin type.

It's versatile, with single and double-pulse modes at 532nm and 1064nm wavelengths, respectively, and adjustable energy settings. The spot size can be set as small as 2mm, or as large as 10mm.

At max power, the MediYag's 1300 mJ of effective energy is enough to complete a full facial skin rejuvenation in under a minute.

A Q-Switched Nd:YAG Without The Gimmicks.

To be perfectly honest, the MediYag is not a revolutionary device. It has no single feature that's going to completely upend the aesthetic world.

But what is there has been painstakingly polished to perfection.

We wanted something easier to use, quicker to master. Something comfortable to operate, day in and day out. Something that lasts.

The MediYag is our ideal vision of what a Q-Switched Nd:YAG laser should be.



It's powerful. Advanced. And it never has any arbitrary software locks or requires branded consumables.

Ask us about the MediYag, and discover one of the industry's most refined answers to tattoo removal, skin rejuvenation, and hyperpigmentation issues.

Every Medicreations Device Is Backed By An Industry-Leading 2-Year “We’ll Fix It” Warranty.

In the rare event your devices malfunction or break due to a technical defect—we’ll **take care of it.**

The benefit of having our manufacturing facility in the U.S. is that we always have access to the very same engineers and technicians who designed and built your machine in the first place. We know exactly how to get your devices working perfectly again. And if we don’t? We’ll replace it at no charge. **Simple.**

We stand behind our products for a full two years—double that of the industry standard. Nor do we require you to pay tens of thousands of dollars for that warranty.

If you buy any Medicreations device, you just won’t have to worry about it for a full 730 days, with options past that.



The Whole Is Greater Than The Sum Of Its Parts, And **Definitely** One Part Alone.

To most med spa operators, a single device can seem like an incredible investment. Here's the calculus.

If a machine works out well, the buyer might earn a return in 2-3 years. At which point they might buy a different one. Rinse and repeat, and they could one day provide the full range of experiences their patients are looking for. Theoretically.

Let's skip that.

Medicreations offers a whole suite of products to deliver comprehensive aesthetic treatments for patients, all at the cost of a single device from other manufacturers. You can provide hair removal alongside skin tightening, body sculpting with muscle toning, or any other synergistic combination. No waiting for years—you could own them all tomorrow.

The difference between 1 device & several is incalculable.

It's a single instrument against an orchestra, a solo competing with a symphony. Only when all the pieces meld together will you find yourself producing that cohesive experience that is "just so."

We'll Teach You Everything About Your Device (And More Besides).

We *love* helping people learn the science and use behind our products. Seriously. So much so we actually built a little academy here in Vegas to do it better.

You'll receive a comprehensive onboarding process when you get any of our lineup at no additional cost. Everything you need to know about usage and care of your devices will be included. That's expected, a given. But we offer more.

**Regular Webinars With Top
Doctors On The Science**

**In-Person Seminars
With Use-Cases**

**Live
Demonstrations**

**Q&A
Sessions**

**Marketing
Help**

Whatever you want to learn, we'll provide.

See If A Medicreations Device Is Perfect For You

Lets talk. We'll make sure you understand the science behind the device you're interested in, who it's ideal for, and how to get the most from it. And if during the process, you find out you're not so sure it's for you?

That's OK.

We only want to sell devices that we sincerely believe will be a benefit to your practice, help you hit the goals you want to hit, and blow your ROI out of the water.

Something that checks every box for you in particular, with NO exceptions.

Maybe that something will be a Medicreations device.

Contact:

Private Practice Solutions LLC

www.ppslasers.com

E: info@ppslasers.com | T: (323) 389-6028

Disclaimer: No representation or warranty will be implied from any description of or claims regarding the equipment or its effectiveness or the ability to achieve any particular clinical results, whether written or oral, as contained in this brochure or similar statements made or furnished to buyer by any person including Medicreations personnel. The use of the equipment requires the expertise of sound professional judgment and the treatments results may be based on operator skill and experience, patient/client suitability, patient/client response to treatment and other factors beyond the control of Medicreations. Operation of the equipment must be in accordance with all applicable laws and regulations and medical and treatments guidelines and for ensuring that each operator of the equipment is adequately trained and qualified to use and operate the equipment safely and properly and to perform procedures in accordance with such law, regulations and guidelines. Medicreations makes no representation or warranties regarding federal, state or local laws or regulations, or medical or treatments guidelines that may apply to the use and operation of the equipment.