

MC MEDICREATIONS

# The MediCO



**Greater Flexibility. Less Upkeep. The Next Stage In  
Fractional CO2 Facial Resurfacing Technology.**



## The Premier Aesthetic Tool To Fight Serious Wrinkles, Sun Damage, Or Discolored Skin.

Sometimes, lightweight treatments don't cut it.

A more intense solution is needed for thick scar tissue, heavily sun-damaged skin, or wrinkles that have had years to "set in."

That's when you reach for the MediCO.

Like other CO2 lasers, the MediCO works by using light to apply carefully controlled damage to the skin, stimulating the body's natural regenerative processes to grow new, healthy, layers.

This gives the patient as close to a "facial reset" as is, today, possible.

Unlike other CO2 lasers, it has a few advantages.

The standard glass tubes for holding the gas are replaced by RF-excited steel canisters, stabilizing energy output and providing a multi-fold increase in durability. Air cooling means the device is whisper-quiet and efficient when in use. And in terms of base specifications, the MediCO boasts the power and flexibility needed to perform even the trickiest treatments.

# America's Leading Clinics Use Medicreation Devices.



## The MediCO:

<b>Touch Screen</b>	10" touch screen LCD monitor
<b>Laser Tube</b>	RF Excited Metal Tube
<b>Scan Depth</b>	Up to 3,000 microns
<b>Fractionated</b>	Yes

# Sick of the Status Quo? So Are We.

The Medicreations team has nearly

# 50

years of combined experience immersed in some of the world's leading aesthetic device brands.



We've learned how things are "supposed to go" from the giants of the industry. **It's terrible.**



**Devices built around...** Expensive, single-use consumables. Manipulative sales tactics. Absurd pricing structures designed to gouge every last cent from care providers.



**So...** We got together some pretty smart scientists, engineers, and business people to do things differently.

## **Our revolutionary idea:**

Make a living *with* the people who use our devices, not *from* them.

# Quality Devices + Transparent Process + Fair Pricing = Our Promise

Aesthetic care professionals should have access to reliable, high-performance devices at prices that actually make sense.

That's why we design our machines to utilize *as few consumables* as possible, even when competitors don't. Why we have a *industry-leading two-year guarantee* for all parts and damage. Why even though we'd put our devices against nearly any other brand's in terms of performance, *we don't charge a premium*.

As for our process, our representatives will never sell you a device that you do not need. If it won't improve your business and make your life richer, it doesn't go out the door. Period.

**“The way it's always been”** isn't the right way. We want to do better.

# Always Built In Our FDA-Compliant Facility Right Here In Las Vegas, NV.

Medcreations devices are some of the only aesthetic care devices fully manufactured in the United States, with at least

# 61%

of all components sourced from U.S. suppliers.

Our FDA-compliant facility features a clean room and full-service center, staffed by our excellent team of professionals. Every device is meticulously built and subject to strenuous quality control checks before shipping out.

We film videos of our process and offer guided tours of our facility on request. We are proud of how we build our devices and what goes into them.

If you're ever down in Vegas, let us know and we'll show you around!

# The Top Characteristic Practitioners Should Look For In A Facial Resurfacing Laser.

Despite what you may hear, it's not simply power.

True, a laser needs enough "weight" to get the job done.

And the MediCO has nothing to be shy about in that regard. In fact, it's 50% stronger than your standard CO2 laser at 30 W, with an adjustable spot size up to 20 x 20. Most are just 12 x 12.

But what truly makes it the ideal practitioner's tool is the degree of customization it allows in each treatment.

Flexibility and judgment, most experienced professionals would agree, is the key to minimizing treatment risk while maximizing efficacy. Analyzing a particular patient's situation and changing the treatment to be more aggressive (or less).

With the MediCO, you can choose from ablative to non-ablative settings. Modify how far the ablative channels are from one another. Adjust the energy output. Switch scan patterns.

This allows a completely customized treatment for each and every patient, ensuring the highest chance for amazing results.

# An Evolution Of The Traditional CO<sub>2</sub>-Powered Laser.

CO<sub>2</sub> lasers are powered by gas, usually stored in a glass tube. Electricity is run through that gas, producing the light necessary for treatment.

But there's a problem with the traditional setup.

This method is relatively crude, can cause instability in light output, and a host of other issues. Among those issues: they need to be re-gassed every 2-3 years of operation assuming moderate usage. A single re-gassing can cost practitioners, on average, \$9,000.

Not a great deal for anyone having to do it.

Instead, the MediCO utilizes steel canisters to store gas. RF waves are used to excite the gas, rather than a direct current of electricity.

This results in light output many times more stable than electricity-excited devices, with fewer instabilities. They last an average of 10-12 years before needed re-gassing – about four times as long.

**How's that for a step forward?**

## Streamlined Design With Air-Cooling And No Counter-Weights.

Standard CO<sub>2</sub> laser designs often include liquid cooling apparatuses. In terms of cooling efficiency, some are quite excellent.

Unfortunately, those tend to be expensive to refill, and introduce additional complexity to the device, which means one more potential point of failure.

The MediCo streamlines the traditional design by incorporating an air-cooling jet instead. It's quieter, and doesn't require the same costly upkeep.

The articulating arm has been specially designed not to require a counter-weight, and moves exactly where you want it to move without any additional complications.

When you need an efficient, potent, flexible solution to facial resurfacing (without the ruinous upkeep costs), consider a MediCO.

# Every Medicreations Device Is Backed By An Industry-Leading 2-Year “We’ll Fix It” Warranty.

In the rare event your devices malfunction or break due to a technical defect—we’ll **take care of it.**

The benefit of having our manufacturing facility in the U.S. is that we always have access to the very same engineers and technicians who designed and built your machine in the first place. We know exactly how to get your devices working perfectly again. And if we don’t? We’ll replace it at no charge. **Simple.**

We stand behind our products for a full two years—double that of the industry standard. Nor do we require you to pay tens of thousands of dollars for that warranty.

If you buy any Medicreations device, you just won’t have to worry about it for a full 730 days, with options past that.

**Contact:**

Private Practice Solutions LLC

[www.ppslasers.com](http://www.ppslasers.com)

E: [info@ppslasers.com](mailto:info@ppslasers.com)

T: (323) 389-6028