

The MediDiode



**3X Faster Laser
Hair Removal.**

**5 Skin
Types.**

**A Fraction Of
The Price.**



At First Glance, The MediDiode Sits Firmly In The “Too Good To Be True” Category.

Cutting-edge 10hz diode laser technology ensures faster, cleaner, less painful hair removal for a wider range of skin types. It's endorsed by renowned professionals in the field of aesthetic care. And an industry-leading two-year device service guarantee covers practically any problem that might unexpectedly ruin your month.

It's half the cost of comparable devices.

When it comes to a premium hair removal investment for your practice that will pay for itself many times over, generate rave reviews from patients, and last for a long, long time—you'd be hard-pressed to find a value that comes close.


America's Leading Clinics Use Medicreation Devices.

"The Medi-Diode is fast, [has a] huge spot size, and allows our patients to get in and out the door in 20 minutes or less with no pain. I would highly recommend this device."

Zach Samios, Owner,
Beverly Hills Rejuvenation Center

The MediDiode:

Skin Types	1-5
Laser Type	Diode
Power	1200 W
Spot Size	12 X 12
Pulse Width	10-800 ms
Frequency	1-10 Hz
Cooling Temp	0-20 C
Avg. Session Length	20 Min
Additional Costs	None
Warranty Period	2 Years
Wavelength	808 nm
Made In	USA



Sick of the Status Quo? So Are We.

The Medicreations team has nearly

50

years of combined experience immersed in some of the world's leading aesthetic device brands.



We've learned how things are "supposed to go" from the giants of the industry. **It's terrible.**



Devices built around... Expensive, single-use consumables. Manipulative sales tactics. Absurd pricing structures designed to gouge every last cent from care providers.



So... We got together some pretty smart scientists, engineers, and business people to do things differently.

Our revolutionary idea:

Make a living *with* the people who use our devices, not *from* them.

Quality Devices + Transparent Process + Fair Pricing = Our Promise

Aesthetic care professionals should have access to reliable, high-performance devices at prices that actually make sense.

That's why we design our machines to utilize *as few consumables* as possible, even when competitors don't. Why we have a *industry-leading two-year guarantee* for all parts and damage. Why even though we'd put our devices against nearly any other brand's in terms of performance, *we don't charge a premium.*

As for our process, our representatives will never sell you a device that you do not need. If it won't improve your business and make your life richer, it doesn't go out the door. Period.

“The way it's always been” isn't the *right* way. We want to do better.

Always Built In Our FDA-Compliant Facility Right Here In Las Vegas, NV.

Medcreations devices are some of the *only* aesthetic care devices fully manufactured in the United States, with at least

61%

of all components sourced from U.S. suppliers.

Our FDA-compliant facility features a clean room and full-service center, staffed by our excellent team of professionals. Every device is meticulously built and subject to strenuous quality control checks before shipping out.


We film videos of our process and offer guided tours of our facility on request. We are proud of how we build our devices and what goes into them.

If you're ever down in Vegas, let us know and we'll show you around!

True High Energy Hair Removal. Faster. Cleaner. Better.

Power matters. When you hook up a 1,200 watt power supply to a 10hz laser tuned to a specific frequency, nearly four times the strength of standard devices, you end up with one of the most effective hair removal tools in existence.

The MediDiode emits energy in a 12x12 mm area, quickly heating up any hair follicle caught inside that square to 40° C in seconds. This saturation essentially “kills” the follicle, making it near-impossible for any new hair to grow. A chilled tip minimizes discomfort throughout the process.



A few sessions with a MediDiode will lead to permanent hair reduction for most patients. And it's thorough.

No More Stamping.

The usual way hair removal is done is that a practitioner activates the laser on a single spot of a person's body, moving in a grid pattern within the treatment area.

Stamp. Move on. Stamp.

This process is tedious, painful for the patient, and requires meticulous attention to detail to avoid missed spots. Skin irritation and the smell of burned hair are standard. For normal lasers this is the only solution – they simply don't have the power to eliminate hair in a quick pass.

Our MediDiode evolves the traditional method. Usage is instead based on rapidly moving the hand device over the target areas, like painting.

The 10hz emitter has enough power that it doesn't need to target a tiny spot to guarantee everything it touches is treated. An XL version with a larger spot size is also available.



That means every treatment can be finished in under 20 minutes, with less effort.

Skin Types 1-5 With One Device.

You shouldn't need two hair removal devices. The MediDiode is a comprehensive solution that can be used for effective hair removal treatment for 5 out of 6 skin types, and it'll work just as well on the darker skin tones as it does on the lighter ones.

That's nearly double the range of comparable devices, and the key is the unique 808nm wavelength utilized by the emitter.

With a MediDiode, right out of the gate, you have a hair removal device that does its job **better, faster, for more people.**

No branded consumables – the only cost is for the machine itself. If you're looking for an all-in-one hair removal solution for a ridiculously low price, a MediDiode might be it.



Type 1

Type 2

Type 3

Type 4

Type 5

Light, Pale White

White, Fair

Medium, White to Olive

Olive, Moderate Brown

Brown, Dark brown

How Much Can You Earn With A MediDiode?

Small Area Treatment- **\$120**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$1200	120	\$14,400
20	\$2400	240	\$28,800
30	\$3600	360	\$43,200
40	\$4800	480	\$57,600
50	\$6000	600	\$72,000

Neck, Upper Lip, Chin, Bikini Line, Armpit, Belly Button, Belly, Back Of Neck, Shoulders, Beard Lines, Cheeks, Hands, Feet, Toes

Large Area Treatment- **\$250**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$2500	120	\$30,000
20	\$5000	240	\$60,000
30	\$7500	360	\$90,000
40	\$10,000	480	\$120,000
50	\$12,500	600	\$150,000

Legs, Chest, Back

Medium Area Treatment- **\$180**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$1800	120	\$21,600
20	\$3600	240	\$43,200
30	\$5400	360	\$64,800
40	\$7200	480	\$86,400
50	\$9000	600	\$108,000

Arms, Bikini, Upper Leg, Lower Leg, Bottom

Full Body Treatment- **\$450**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$4500	120	\$54,000
20	\$9000	240	\$108,000
30	\$13,500	360	\$162,000
40	\$18,000	480	\$216,000
50	\$22,500	600	\$270,000

Every Medicreations Device Is Backed By An Industry-Leading 2-Year “We’ll Fix It” Warranty.

In the rare event your devices malfunction or break due to a technical defect—we’ll **take care of it.**

The benefit of having our manufacturing facility in the U.S. is that we always have access to the very same engineers and technicians who designed and built your machine in the first place. We know exactly how to get your devices working perfectly again. And if we don’t? We’ll replace it at no charge. **Simple.**

We stand behind our products for a full two years—double that of the industry standard. Nor do we require you to pay tens of thousands of dollars for that warranty.

If you buy any Medicreations device, you just won’t have to worry about it for a full 730 days, with options past that.



The Whole Is Greater Than The Sum Of Its Parts, And **Definitely** One Part Alone.

To most med spa operators, a single device can seem like an incredible investment. Here's the calculus.

If a machine works out well, the buyer might earn a return in 2-3 years. At which point they might buy a different one. Rinse and repeat, and they could one day provide the full range of experiences their patients are looking for. Theoretically.

Let's skip that.

Medicreations offers a whole suite of products to deliver comprehensive aesthetic treatments for patients, all at the cost of a single device from other manufacturers. You can provide hair removal alongside skin tightening, body sculpting with muscle toning, or any other synergistic combination. No waiting for years—you could own them all tomorrow.

The difference between 1 device & several is incalculable.

It's a single instrument against an orchestra, a solo competing with a symphony. Only when all the pieces meld together will you find yourself producing that cohesive experience that is "just so."

We'll Teach You Everything About Your Device (And More Besides).

We *love* helping people learn the science and use behind our products. Seriously. So much so we actually built a little academy here in Vegas to do it better.

You'll receive a comprehensive onboarding process when you get any of our lineup at no additional cost. Everything you need to know about usage and care of your devices will be included. That's expected, a given. But we offer more.

Regular Webinars With Top Doctors On The Science

In-Person Seminars With Use-Cases

Live Demonstrations

Q&A Sessions

Marketing Help

Whatever you want to learn, we'll provide.



See If A Medicreations Device Is Perfect For You.

Let talk. We'll make sure you understand the science behind the device you're interested in, who it's ideal for, and how to get the most from it. And if during the process, you find out you're not so sure it's for you?

That's OK.

We only want to sell devices that we sincerely believe will be a benefit to your practice, help you hit the goals you want to hit, and blow your ROI out of the water. Something that checks every box for you in particular, with NO exceptions.

Maybe that something will be a Medicreations device.

Contact:

Private Practice Solutions LLC

www.ppslasers.com

E: info@ppslasers.com | T: (323) 389-6028

Disclaimer: Claims of 3X faster hair removal based on standard 3hz diode laser devices. No representation or warranty will be implied from any description of or claims regarding the equipment or its effectiveness or the ability to achieve any particular clinical results, whether written or oral, as contained in this brochure or similar statements made or furnished to buyer by any person including Medicreations personnel. The use of the equipment requires the expertise of sound professional judgment and the treatments results may be based on operator skill and experience, patient/client suitability, patient/client response to treatment and other factors beyond the control of Medicreations. Operation of the equipment must be in accordance with all applicable laws and regulations and medical and treatments guidelines and for ensuring that each operator of the equipment is adequately trained and qualified to use and operate the equipment safely and properly and to perform procedures in accordance with such law, regulations and guidelines. Medicreations makes no representation or warranties regarding federal, state or local laws or regulations, or medical or treatments guidelines that may apply to the use and operation of the equipment. Published on March 13th, 2020.