

MC MEDICREATIONS

# The MediLight



**Photon Pulse Train™ Technology  
Enables Complete Treatment In One Pass.**



## A True Jack-Of-All-Trades (With A Greater Punch).

When it comes to versatility, the IPL (Intense Pulsed Light) is king.

It treats hyperpigmentation issues.

Uneven tone. Sunspots. Birthmarks.

It can help with acne. Professionals reach for an IPL to smooth out skin texture. Even hair removal is part of its repertoire.

But not all IPLs are created equal.

You need something potent enough to perform treatments quickly, with high efficacy and ease-of-use. Equally important is a priority on minimizing risk and pain for patients.

In that regard, at first glance, the MediLight IPL doesn't seem too different than the competition.

Yet it features unique Photon Pulse Train™ technology that allows for a greater concentration of energy on the intended target, with lower risk of burning the surrounding skin. That means practitioners can perform more aggressive treatment options with greater ease-of-mind. One pass of energy, instead of several.

It is unparalleled when it comes to minimizing upkeep. A special swappable filter system means it'll operate at peak efficiency for many years without service, and the filters won't "burn out" in a few months.

Best of all? No arbitrary locks on modalities. What you see is what you get.

# America's Leading Clinics Use Medicreation Devices.



## The MediLight:

	IPL Mode	LHR Mode
<b>LCD Display</b>	15" color touch LCD screen	
<b>Pulse Train</b>	Yes-up to 6 sub-pulses.	NA
<b>Hertz</b>	3Hz	1-10Hz
<b>Chill Tip</b>	Adjustable 0-20C	0-20C
<b>Energy</b>	1-15 J/cm2	1-15 J/cm2
<b>Pulse Duration</b>	1-9.9 ms	1-50 ms
<b>Spot Size</b>	12x50 mm	12x50 mm
<b>Wavelength Filters</b>	480/530/590/640/690nm	640nm

# Sick of the Status Quo? So Are We.

The Medicreations team has nearly

# 50

years of combined experience immersed in some of the world's leading aesthetic device brands.



We've learned how things are "supposed to go" from the giants of the industry. **It's terrible.**



**Devices built around...** Expensive, single-use consumables. Manipulative sales tactics. Absurd pricing structures designed to gouge every last cent from care providers.



**So...** We got together some pretty smart scientists, engineers, and business people to do things differently.

## **Our revolutionary idea:**

Make a living *with* the people who use our devices, not *from* them.

# Quality Devices + Transparent Process + Fair Pricing = Our Promise

Aesthetic care professionals should have access to reliable, high-performance devices at prices that actually make sense.

That's why we design our machines to utilize ***as few consumables*** as possible, even when competitors don't. Why we have a ***industry-leading two-year guarantee*** for all parts and damage. Why even though we'd put our devices against nearly any other brand's in terms of performance, ***we don't charge a premium.***

As for our process, our representatives will never sell you a device that you do not need. If it won't improve your business and make your life richer, it doesn't go out the door. Period.

**“The way it's always been”** isn't the right way. We want to do better.

# Always Built In Our FDA-Compliant Facility Right Here In Las Vegas, NV.

Medicreations devices are some of the only aesthetic care devices fully manufactured in the United States, with at least

# 61%

of all components sourced from U.S. suppliers.

Our FDA-compliant facility features a clean room and full-service center, staffed by our excellent team of professionals. Every device is meticulously built and subject to strenuous quality control checks before shipping out.

We film videos of our process and offer guided tours of our facility on request. We are proud of how we build our devices and what goes into them.

If you're ever down in Vegas, let us know and we'll show you around!

# Power + Safety. Photon Pulse Train™ Technology Enables Rapid Treatment And Consistent Patient Comfort.


It might surprise you to learn that most IPLs never utilize 100% of their effective energy for any skin treatment.

That's because as the IPL is utilized on things like hair, blood vessels, or pigmentation, the temperature of the skin climbs along with the overall energy used.

Try to do too aggressive a treatment? You risk producing burns.

**Switching the MediLight to Pulse Train mode completely upends this dynamic.** This setting breaks up larger pulses into sub-pulses. As the energy builds, the device takes more frequent “rest” intervals between shorter pulses.

## How Photon Pulse Train™ Works



This keeps the skin cool and comfortable while treatment efficacy on the target hits its peak. It's a no-compromise solution to IPL treatment.

Most MediLight treatments require only a single pass over the face or body, whereas other devices might require several.

# Customize Individual IPL Treatments To An Absurd Degree.

You can modify practically every setting there is for the MediLight. Each one is independent of any other.

You can adjust the duration of each pulse. The pause between them. Even the length of energy emission for each particular sub-pulse.

There are, of course, recommended settings for each and every treatment you'll want to perform. Settings that have been battle-tested in the real world and proven to be effective at what they should do.

But if a practitioner really knows what they're doing, the MediLight's inherent flexibility means they can create completely customized, one-of-a-kind treatments for their patients.

## “Swap-In” Filter System Means Greater Reliability, Less Burned-Out Filters.

The bane of many IPL device owners is just how quickly the filters burn out after repeated treatments.

For some, this can take years. For the truly unlucky – months.

And it can take thousands (even tens of thousands) of dollars to replace these filters.

We've gotten rid of that gamble with a completely different system for the MediLight. It does come with one pre-coated sapphire handpiece for hair removal.

But the other handpiece is an uncoated IPL emitter, with a swappable filter system. These filters are located further away from the bulbs than a “traditional” system. This significantly reduces the risk of anything burning out, and allows more reliable and longer-lasting usage.

**And not only that.**

# No Additional Add-Ons. One Handpiece With Every Filter You'll Ever Need, Included.

Most popular IPL devices come with a few pre-enabled modalities.


Then, one by one, you buy the individual pre-coated handpieces you need to unlock different treatment options. Eventually, you can actually perform all the treatments you want to (after spending a small fortune).

We have a different approach.

**We just never lock anything in the first place.**

Every MediLight IPL comes equipped with *everything* you'll need to do any treatment it can do. Immediately. With no add-ons or fees.

## How The MediLight Filters Work:



Simply use one of the five swap-in filters or the pre-coated hair removal handpiece and you're good to go.

When you need versatility, reliability, and greater efficiency of power in an IPL: a MediLight might be what you're looking for.

When you need versatility, reliability, and greater efficiency of power in an IPL: **a MediLight might be what you're looking for.**

# How Much Can You Earn With A MediLight?

Small Area Treatment — **\$120 /**

Package — **\$950 (10 sessions)**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$1200	120	\$14,400
20	\$2400	240	\$28,800
30	\$3600	360	\$43,200
40	\$4800	480	\$57,600
50	\$6000	600	\$72,000

Neck, Upper Lip, Chin, Bikini Line, Armpit, Belly button, Belly, Back of Neck, Shoulders, Beard lines, Cheeks, Hands, Feet, Toes

Large Area Treatment — **\$250 /**

Package — **\$1950 (10 sessions)**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$2500	120	\$30,000
20	\$5000	240	\$60,000
30	\$7500	360	\$90,000
40	\$10,000	480	\$120,000
50	<b>\$12,500</b>	600	<b>\$150,000</b>

Legs, Chest, Back

Medium Area Treatment — **\$180 /**

Package — **\$1450 (10 sessions)**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$1800	120	\$21,600
20	\$3600	240	\$43,200
30	\$5400	360	\$64,800
40	\$7200	480	\$86,400
50	<b>\$9000</b>	600	<b>\$108,000</b>

Arms, Bikini, Upper Leg, Lower Leg, Bottom

Full Body — **\$450 /**

Package — **\$3700 (10 sessions)**

Monthly Treatments	Monthly Revenue	Annual Treatments	Annual Revenue
10	\$4500	120	\$54,000
20	\$9000	240	\$108,000
30	\$13,500	360	\$162,000
40	\$18,000	480	\$216,000
50	\$22,500	600	\$270,000

# Every Medicreations Device Is Backed By An Industry-Leading 2-Year “We’ll Fix It” Warranty.

In the rare event your devices malfunction or break due to a technical defect—we’ll **take care of it.**

The benefit of having our manufacturing facility in the U.S. is that we always have access to the very same engineers and technicians who designed and built your machine in the first place. We know exactly how to get your devices working perfectly again. And if we don’t? We’ll replace it at no charge. **Simple.**

We stand behind our products for a full two years—double that of the industry standard. Nor do we require you to pay tens of thousands of dollars for that warranty.

If you buy any Medicreations device, you just won’t have to worry about it for a full 730 days, with options past that.



## The Whole Is Greater Than The Sum Of Its Parts, And Definitely One Part Alone.

To most med spa operators, a single device can seem like an incredible investment. Here's the calculus.

If a machine works out well, the buyer might earn a return in 2-3 years. At which point they might buy a different one. Rinse and repeat, and they could one day provide the full range of experiences their patients are looking for. Theoretically.

*Let's skip that.*

Medicreations offers a whole suite of products to deliver comprehensive aesthetic treatments for patients, all at the cost of a single device from other manufacturers. You can provide hair removal alongside skin tightening, body sculpting with muscle toning, or any other synergistic combination. No waiting for years—you could own them all tomorrow.

**The difference between 1 device & several is incalculable.**

It's a single instrument against an orchestra, a solo competing with a symphony. Only when all the pieces meld together will you find yourself producing that cohesive experience that is "just so."

# We'll Teach You Everything About Your Device (And More Besides).

We *love* helping people learn the science and use behind our products. Seriously. So much so we actually built a little academy here in Vegas to do it better.

You'll receive a comprehensive onboarding process when you get any of our lineup at no additional cost. Everything you need to know about usage and care of your devices will be included. That's expected, a given. But we offer more.

**Regular Webinars With Top Doctors On The Science**

**In-Person Seminars With Use-Cases**

**Live Demonstrations**

**Q&A Sessions**

**Marketing Help**

Whatever you want to learn, we'll provide.

# See If A Medicreations Device Is Perfect For You

**Lets talk.** We'll make sure you understand the science behind the device you're interested in, who it's ideal for, and how to get the most from it. And if during the process, you find out you're not so sure it's for you?

**That's OK.**

We only want to sell devices that we sincerely believe will be a benefit to your practice, help you hit the goals you want to hit, and blow your ROI out of the water.

Something that checks every box for you in particular, with NO exceptions.

**Maybe that something will be a Medicreations device.**

## Contact:

Private Practice Solutions LLC

[www.ppslasers.com](http://www.ppslasers.com)

E: [info@ppslasers.com](mailto:info@ppslasers.com) | T: (323) 389-6028

**Disclaimer:** No representation or warranty will be implied from any description of or claims regarding the equipment or its effectiveness or the ability to achieve any particular clinical results, whether written or oral, as contained in this brochure or similar statements made or furnished to buyer by any person including Medicreations personnel. The use of the equipment requires the expertise of sound professional judgment and the treatments results may be based on operator skill and experience, patient/client suitability, patient/client response to treatment and other factors beyond the control of Medicreations. Operation of the equipment must be in accordance with all applicable laws and regulations and medical and treatments guidelines and for ensuring that each operator of the equipment is adequately trained and qualified to use and operate the equipment safely and properly and to perform procedures in accordance with such law, regulations and guidelines. Medicreations makes no representation or warranties regarding federal, state or local laws or regulations, or medical or treatments guidelines that may apply to the use and operation of the equipment.